PUL WASHABLE UNDERPAD

A PUL(PU Laminating) Washable Underpad(Bed Pad & Chair Pad) is a piece

of soft cloth that is quilted with a absorbent nonwoven underneath it and is laminated to a waterproof PU film, and then laminated to a non-slip backing again. The purpose of the product is to capture urine and absorb it into the nonwoven material so it does not ruin furniture. A PUL washable underpad is better for the environment than a PVC washable underpad, and it should have 4 layers at least: Surface, Soaker, Barrier and Backing.

1. The normal materials for the 4 layers are as below:

Surface: TWILL, Poplin, Jersey, Terry, 100%cotton Flannel, 100% brushed poly Knit, 100%poly Microfibre, etc..

Soaker: Poly/Rayon blend nonwoven, 6 ounce, 8 ounce, 10 ounce, 12 ounce, etc...

Barrier: PU film, 0.02mm-0.05mm Backing: 100%poly knit, non-slip

2. Absorbency level suggestions:

Mild incontinence: 6 ounce or less soaker

Moderate incontinence: 8 ounce or 10 ounce soaker

Heavy incontinence: 12 ounce or more soaker

3. Features:

- 4-layer design offers perfect protection for your skin and bedding
- · Quilted construction ensures many times washing and maintains its shape
- · Washable over 70 times at 90°C water
- Tuck-in flaps are available as requested
- Non-slip bottom layer keeps the pad in place
- · Better for the environment than disposables

























A PVC Washable Underpad(Bed Pad & Chair Pad) is a piece of soft cloth that is quilted with a absorbent nonwoven underneath it and is sewn to a waterproof vinyl backing. The purpose of the product is to capture urine and absorb it into the nonwoven material so it does not ruin furniture. So a standard PVC washable underpad should have 3 layers at least: Surface, Soaker and Barrier.

1. The normal materials for the 3 layers are as below:

Surface: IBEX, Birdeyes, TWILL, Poplin, Jersey, Terry, 100%cotton Flannel, 100% brushed poly Knit, 100%poly

Microfibre, etc..

Soaker: Poly/Rayon blend nonwoven, 6 ounce, 8 ounce, 10 ounce, 12 ounce, etc...

Barrier: 100%vinyl extruded to 100%poly knit, non-slip backing

2. Absorbency level suggestions:

Mild incontinence: 6 ounce or less soaker


Moderate incontinence: 8 ounce or 10 ounce soaker


Heavy incontinence: 12 ounce or more soaker

3. Features:

- 3-layer design offers perfect protection for your skin and bedding
- Quilted construction ensures many times washing and maintains its shape
- · Medically grade vinvl is harmless to your health
- · Washable over 150 times at 90°C water
- Tuck-in flaps are available as requested
- Non-slip bottom layer keeps the pad in place
- Better for the environment than disposables









PBP-109 Size: 80X90cm other customized sizes.



PBP-110 Size: 80X90cm other customized sizes.



PBP-111 Size: 80X90cm other customized sizes.



